




ACTRT2 Antibody

Product Code	CSB-PA001256LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDY3
Immunogen	Recombinant Human Actin-related protein T2 protein (1-376AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein T2 (ARP-T2) (Actin-related protein M2), ACTRT2, ARPM2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACTRT2

Image



Western blot

All lanes: ACTRT2 antibody at 6µg/ml

Lane 1: HepG2 whole cell lysate

Lane 2: A431 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 42 kDa

Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.