



ADAM10 Antibody, Biotin conjugated

Product Code	CSB-PA001273LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14672
Immunogen	Recombinant Human Disintegrin and metalloproteinase domain-containing protein 10 protein (214-504AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disintegrin and metalloproteinase domain-containing protein 10 (ADAM 10) (EC 3.4.24.81) (CDw156) (Kuzbanian protein homolog) (Mammalian disintegrin-metalloprotease) (CD antigen CD156c), ADAM10, KUZ MADM
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ADAM10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.