



ADAMTS7 Antibody, FITC conjugated

Product Code	CSB-PA891951LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKP4
Immunogen	Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 7 protein (242-593AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Metalloprotease that may play a role in the degradation of COMP.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	A disintegrin and metalloproteinase with thrombospondin motifs 7 (ADAM-TS 7) (ADAM-TS7) (ADAMTS-7) (EC 3.4.24.-) (COMPase), ADAMTS7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADAMTS7