



ADAP2 Antibody

Product Code	CSB-PA889071LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPF8
Immunogen	Recombinant Human Arf-GAP with dual PH domain-containing protein 2 protein (110-381AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arf-GAP with dual PH domain-containing protein 2 (Centaurin-alpha-2) (Cnt-a2), ADAP2, CENTA2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADAP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.