



ADGRL4 Antibody, Biotin conjugated

Product Code	CSB-PA872541LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBW9
Immunogen	Recombinant Human Adhesion G protein-coupled receptor L4 protein (85-299AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adhesion G protein-coupled receptor L4 (EGF, latrophilin and seven transmembrane domain-containing protein 1) (EGF-TM7-latrophilin-related protein) (ETL protein), ADGRL4, ELTD1 ETL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADGRL4