



ADGRL4 Antibody, Biotin conjugated

Product Code	CSB-PA872541LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBW9
Immunogen	Recombinant Human Adhesion G protein-coupled receptor L4 protein (85-299AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adhesion G protein-coupled receptor L4 (EGF, latrophilin and seven transmembrane domain-containing protein 1) (EGF-TM7-latrophilin-related protein) (ETL protein), ADGRL4, ELTD1 ETL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADGRL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.