



ADPRHL1 Antibody

Product Code	CSB-PA839860LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDY3
Immunogen	Recombinant Human [Protein ADP-ribosylarginine] hydrolase-like protein 1 protein (1-272AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	[Protein ADP-ribosylarginine] hydrolase-like protein 1 (EC 3.2.-.-) (ADP-ribosylhydrolase 2), ADPRHL1, ARH2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ADPRHL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.