





AGBL2 Antibody

Product Code	CSB-PA732978LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5U5Z8
Immunogen	Recombinant Human Cytosolic carboxypeptidase 2 protein (10-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	AGBL 2 antibody; AGBL2 antibody; ATP/GTP binding protein like 2 antibody; ATP/GTP-binding protein-like 2 antibody; Carboxypeptidase 2 cytosolic antibody; CBPC2_HUMAN antibody; CCP2 antibody; Cytosolic carboxypeptidase 2 antibody; FLJ23598 antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AGBL2

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA732978LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA732978LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.