



AGBL2 Antibody, Biotin conjugated

Product Code	CSB-PA732978LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5U5Z8
Immunogen	Recombinant Human Cytosolic carboxypeptidase 2 protein (10-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic carboxypeptidase 2 (EC 3.4.17.-) (ATP/GTP-binding protein-like 2), AGBL2, CCP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AGBL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.