



# AGR3 Antibody

<b>Product Code</b>	CSB-PA848404LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8TD06
<b>Immunogen</b>	Recombinant Human Anterior gradient protein 3 protein (22-166AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Anterior gradient protein 3 (AG-3) (AG3) (hAG-3) (Anterior gradient 3 homolog) (Breast cancer membrane protein 11) (Protein disulfide isomerase family A, member 18), AGR3, BCMP11 PDIA18
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cancer
<b>Target Names</b>	AGR3