



AGR3 Antibody

Product Code	CSB-PA848404LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD06
Immunogen	Recombinant Human Anterior gradient protein 3 protein (22-166AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anterior gradient protein 3 (AG-3) (AG3) (hAG-3) (Anterior gradient 3 homolog) (Breast cancer membrane protein 11) (Protein disulfide isomerase family A, member 18), AGR3, BCMP11 PDIA18
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	AGR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.