




ANKRD10 Antibody

Product Code	CSB-PA878915LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXR5
Immunogen	Recombinant Human Ankyrin repeat domain-containing protein 10 protein (1-220AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ankyrin repeat domain-containing protein 10, ANKRD10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ANKRD10

Image



Western Blot

Positive WB detected in: Mouse heart tissue, Mouse kidney tissue, Mouse brain tissue

All lanes: ANKRD10 antibody at 4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45, 24 kDa

Observed band size: 50 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.