



# AP1S1 Antibody

<b>Product Code</b>	CSB-PA001866LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P61966
<b>Immunogen</b>	Recombinant Human AP-1 complex subunit sigma-1A protein (1-127AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	AP-1 complex subunit sigma-1A (Adaptor protein complex AP-1 subunit sigma-1A) (Adaptor-related protein complex 1 subunit sigma-1A) (Clathrin assembly protein complex 1 sigma-1A small chain) (Clathrin coat assembly protein AP19) (Golgi adaptor HA1/AP1 adaptin sigma-1A subunit) (HA1 19 kDa subunit) (Sigma 1a subunit of AP-1 clathrin) (Sigma-adaptin 1A) (Sigma1A-adaptin), AP1S1, AP19 CLAPS1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	AP1S1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.