




AP3S2 Antibody

Product Code	CSB-PA001880LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59780
Immunogen	Recombinant Human AP-3 complex subunit sigma-2 protein (1-193AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	AP-3 complex subunit sigma-2 (AP-3 complex subunit sigma-3B) (Adaptor-related protein complex 3 subunit sigma-2) (Sigma-3B-adaptin) (Sigma3B-adaptin) (Sigma-adaptin 3b), AP3S2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	AP3S2

Image



Western blot
 All lanes: AP3S2 antibody at 6µg/ml + A549 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 23, 44, 15 kDa
 Observed band size: 23 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.