



APOB Antibody

Product Code	CSB-PA001918LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04114
Immunogen	Recombinant Human Apolipoprotein B-100 protein (28-127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apolipoprotein B-100 (Apo B-100) [Cleaved into: Apolipoprotein B-48 (Apo B-48)], APOB
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	APOB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.