





ARHGAP25 Antibody

Product Code	CSB-PA002026LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42331
Immunogen	Recombinant Human Rho GTPase-activating protein 25 protein (17-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho GTPase-activating protein 25 (Rho-type GTPase-activating protein 25), ARHGAP25, KIAA0053
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGAP25

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA002026LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA002026LA01HU at dilution of 1:100