





ARMCX2 Antibody

Product Code	CSB-PA755485LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L311
Immunogen	Recombinant Human Armadillo repeat-containing X-linked protein 2 protein (161-460AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Armadillo repeat-containing X-linked protein 2 (ARM protein lost in epithelial cancers on chromosome X 2) (Protein ALEX2), ARM CX2, ALEX2 KIAA0512
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	ARMCX2


Image



Western blot
 All lanes: ARM CX2 antibody at 5µg/ml + MCF-7 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 66 kDa
 Observed band size: 66 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA755485LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA755485LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.