




# ARRDC3 Antibody

<b>Product Code</b>	CSB-PA839293LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96B67
<b>Immunogen</b>	Recombinant Human Arrestin domain-containing protein 3 protein (1-270AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Arrestin domain-containing protein 3 (TBP-2-like inducible membrane protein) (TLIMP), ARRDC3, KIAA1376
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology

## Image



### Western Blot

Positive WB detected in: MCF-7 whole cell lysate, Mouse skeletal muscle tissue

All lanes: ARRDC3 antibody at 4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 47 kDa

Observed band size: 47 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.