




ASPHD2 Antibody

Product Code	CSB-PA760787LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ICH7
Immunogen	Recombinant Human Aspartate beta-hydroxylase domain-containing protein 2 protein (80-369AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aspartate beta-hydroxylase domain-containing protein 2 (EC 1.14.11.-), ASPHD2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ASPHD2

Image



Western blot

All lanes: ASPHD2 antibody at 4µg/ml

Lane 1: Hela whole cell lysate

Lane 2: HL60 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 42 kDa

Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.