



ATXN7L3 Antibody, HRP conjugated

Product Code	CSB-PA613538LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14CW9
Immunogen	Recombinant Human Ataxin-7-like protein 3 protein (51-347AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ataxin-7-like protein 3 (SAGA-associated factor 11 homolog), ATXN7L3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ATXN7L3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.