




BAIAP2L1 Antibody

Product Code	CSB-PA887045LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHR4
Immunogen	Recombinant Human Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 1 protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 1 (BAI1-associated protein 2-like protein 1) (Insulin receptor tyrosine kinase substrate), BAIAP2L1, IRTKS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BAIAP2L1

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA887045LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.