





BDH2 Antibody

Product Code	CSB-PA874814LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUT1
Immunogen	Recombinant Human 3-hydroxybutyrate dehydrogenase type 2 protein (1-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-hydroxybutyrate dehydrogenase type 2 (EC 1.1.1.-) (EC 1.1.1.30) (Dehydrogenase/reductase SDR family member 6) (Oxidoreductase UCPA) (R-beta-hydroxybutyrate dehydrogenase) (Short chain dehydrogenase/reductase family 15C member 1), BDH2, DHRS6 SDR15C1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	BDH2

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA874814LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA874814LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.