




# BTN2A1 Antibody

<b>Product Code</b>	CSB-PA755482LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q7KYR7
<b>Immunogen</b>	Recombinant Human Butyrophilin subfamily 2 member A1 protein (29-248AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:300, IF:1:50-1:300
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	BTN2A1

## Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA755482LA01HU at dilution of 1:300



Immunofluorescent analysis of U251 cells using CSB-PA755482LA01HU at dilution of 1:300 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)