





BTN2A1 Antibody

Product Code	CSB-PA755482LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7KYR7
Immunogen	Recombinant Human Butyrophilin subfamily 2 member A1 protein (29-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:300, IF:1:50-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BTN2A1

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA755482LA01HU at dilution of 1:300



Immunofluorescent analysis of U251 cells using CSB-PA755482LA01HU at dilution of 1:300 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.