



BTN2A1 Antibody, FITC conjugated

Product Code	CSB-PA755482LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7KYR7
Immunogen	Recombinant Human Butyrophilin subfamily 2 member A1 protein (29-248AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Butyrophilin subfamily 2 member A1, BTN2A1, BT2.1 BTF1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BTN2A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.