



BTN3A2 Antibody

Product Code	CSB-PA002874LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78410
Immunogen	Recombinant Human Butyrophilin subfamily 3 member A2 protein (30-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Butyrophilin subfamily 3 member A2, BTN3A2, BT3.2 BTF3 BTF4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	BTN3A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.