



C12orf57 Antibody

Product Code	CSB-PA859105LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99622
Immunogen	Recombinant Human Protein C10 protein (2-126AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein C10, C12orf57, C10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C12orf57
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.