



C12orf57 Antibody

Product Code	CSB-PA859105LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99622
Immunogen	Recombinant Human Protein C10 protein (2-126AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein C10, C12orf57, C10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C12orf57