





C15orf40 Antibody

Product Code	CSB-PA837852LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUR7
Immunogen	Recombinant Human UPF0235 protein C15orf40 protein (28-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UPF0235 protein C15orf40, C15orf40
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C15orf40

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA837852LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA837852LA01HU at dilution of 1:100