



C1orf54 Antibody, FITC conjugated

Product Code	CSB-PA840992LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWF1
Immunogen	Recombinant Human Uncharacterized protein C1orf54 protein (17-131AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C1orf54, C1orf54
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C1orf54
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.