



FNDC11 Antibody, HRP conjugated

Product Code	CSB-PA883611LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BVV2
Immunogen	Recombinant Human Uncharacterized protein C20orf195 protein (10-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibronectin type III domain-containing protein 11, FNDC11, C20orf195
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FNDC11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.