



C7orf33 Antibody

Product Code	CSB-PA840977LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WU49
Immunogen	Recombinant Human Uncharacterized protein C7orf33 protein (1-177AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C7orf33, C7orf33
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C7orf33
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.