



# C7orf33 Antibody

Product Code	CSB-PA840977LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WU49
Immunogen	Recombinant Human Uncharacterized protein C7orf33 protein (1-177AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C7orf33, C7orf33
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C7orf33