




C8orf48 Antibody

Product Code	CSB-PA853456LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LL4
Immunogen	Recombinant Human Uncharacterized protein C8orf48 protein (1-319AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C8orf48, C8orf48
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C8orf48

Image



Western blot

All lanes: C8orf48 antibody at 1µg/ml

Lane 1: HepG2 whole cell lysate

Lane 2: Mouse liver tissue

Lane 3: MCF-7 whole cell lysate


Lane 4: K562 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 37 kDa

Observed band size: 37 kDa



Immunohistochemistry of paraffin-embedded

human placenta tissue using CSB-

PA853456LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.