





CA11 Antibody

Product Code	CSB-PA004366LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75493
Immunogen	Recombinant Human Carbonic anhydrase-related protein 11 protein (24-328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carbonic anhydrase-related protein 11 (CA-RP XI) (CA-XI) (CARP XI) (Carbonic anhydrase-related protein 2) (CA-RP II) (CARP-2), CA11, CARP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CA11

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA004366LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human stomach tissue using CSB-PA004366LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.