



CAB39L Antibody, HRP conjugated

Product Code	CSB-PA872527LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9S4
Immunogen	Recombinant Human Calcium-binding protein 39-like protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding protein 39-like (Antigen MLAA-34) (MO25beta) (Mo25-like protein), CAB39L
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	CAB39L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.