



Camk2n1 Antibody

Product Code	CSB-PA004469HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6QWF9
Immunogen	Recombinant Mouse Calcium/calmodulin-dependent protein kinase II inhibitor 1 protein (1-78AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium/calmodulin-dependent protein kinase II inhibitor 1 (calcium/calmodulin-dependent protein kinase II inhibitor alpha) (mCaMKIINalpha), Camk2n1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Camk2n1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.