



Camk2n1 Antibody

Product Code	CSB-PA004469HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6QWF9
Immunogen	Recombinant Mouse Calcium/calmodulin-dependent protein kinase II inhibitor 1 protein (1-78AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium/calmodulin-dependent protein kinase II inhibitor 1 (calcium/calmodulin-dependent protein kinase II inhibitor alpha) (mCaMKIINalpha), Camk2n1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Camk2n1