



SGF29 Antibody, Biotin conjugated

Product Code	CSB-PA856928LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96ES7
Immunogen	Recombinant Human SAGA-associated factor 29 homolog protein (1-293AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SAGA-associated factor 29 (Coiled-coil domain-containing protein 101) (SAGA complex-associated factor 29), SGF29, CCDC101
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SGF29
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.