



SGF29 Antibody, Biotin conjugated

Product Code	CSB-PA856928LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96ES7
Immunogen	Recombinant Human SAGA-associated factor 29 homolog protein (1-293AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SAGA-associated factor 29 (Coiled-coil domain-containing protein 101) (SAGA complex-associated factor 29), SGF29, CCDC101
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SGF29