





MCUB Antibody

Product Code	CSB-PA873671LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWR8
Immunogen	Recombinant Human Calcium uniporter regulatory subunit MCUB, mitochondrial protein (89-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium uniporter regulatory subunit MCUB, mitochondrial (MCUB) (Coiled-coil domain-containing protein 109B), MCUB, CCDC109B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MCUB

Image



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA873671LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA873671LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)