




CCDC112 Antibody

Product Code	CSB-PA818774LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEF3
Immunogen	Recombinant Human Coiled-coil domain-containing protein 112 protein (201-446AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coiled-coil domain-containing protein 112 (Mutated in bladder cancer protein 1), CCDC112, MBC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCDC112

Image



Western blot

All lanes: CCDC112 antibody at 12µg/ml

Lane 1: Jurkat whole cell lysate

Lane 2: HeLa whole cell lysate

Lane 3: HepG2 whole cell lysate


Lane 4: A549 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 54, 62 kDa


Observed band size: 54 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA818774LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA818774LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA818774LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.