




CCDC127 Antibody

Product Code	CSB-PA842643LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BQ5
Immunogen	Recombinant Human Coiled-coil domain-containing protein 127 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coiled-coil domain-containing protein 127, CCDC127
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCDC127

Image



Western blot
 All lanes: CCDC127 antibody at 6µg/ml + Mouse kidney tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 31 kDa
 Observed band size: 31 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.