





CCDC24 Antibody

Product Code	CSB-PA836665LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N4L8
Immunogen	Recombinant Human Coiled-coil domain-containing protein 24 protein (1-307AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coiled-coil domain-containing protein 24, CCDC24
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CCDC24

Image



Western blot
All lanes: CCDC24 antibody at 6µg/ml + MCF-7 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 35, 31 kDa
Observed band size: 35 kDa



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA836665LA01HU at dilution of 1:100