



WASHC3 Antibody, Biotin conjugated

Product Code	CSB-PA896504LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3C0
Immunogen	Recombinant Human WASH complex subunit 3 protein (1-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WASH complex subunit 3 (Coiled-coil domain-containing protein 53), WASHC3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WASHC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.