



CEACAM1 Antibody, HRP conjugated

Product Code	CSB-PA10749B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13688
Immunogen	Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 1 protein (35-428AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carcinoembryonic antigen-related cell adhesion molecule 1 (Biliary glycoprotein 1) (BGP-1) (CD antigen CD66a), CEACAM1, BGP BGP1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CEACAM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.