



CHID1 Antibody, FITC conjugated

Product Code	CSB-PA883614LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWS9
Immunogen	Recombinant Human Chitinase domain-containing protein 1 protein (20-393AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chitinase domain-containing protein 1 (Stabilin-1-interacting chitinase-like protein) (SI-CLP), CHID1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CHID1