



CHST10 Antibody

Product Code	CSB-PA005403LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43529
Immunogen	Recombinant Human Carbohydrate sulfotransferase 10 protein (28-356AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carbohydrate sulfotransferase 10 (EC 2.8.2.-) (HNK-1 sulfotransferase) (HNK-1ST) (HNK1ST) (HuHNK-1ST), CHST10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CHST10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.