



CLEC1A Antibody, FITC conjugated

Product Code	CSB-PA843319LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NC01
Immunogen	Recombinant Human C-type lectin domain family 1 member A protein (74-280AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-type lectin domain family 1 member A (C-type lectin-like receptor 1) (CLEC-1), CLEC1A, CLEC1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CLEC1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.