





CLEC4C Antibody

Product Code	CSB-PA855470LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WTT0
Immunogen	Recombinant Human C-type lectin domain family 4 member C protein (45-213AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-type lectin domain family 4 member C (Blood dendritic cell antigen 2) (BDCA-2) (C-type lectin superfamily member 7) (Dendritic lectin) (CD antigen CD303), CLEC4C, BDCA2 CLECSF11 CLECSF7 DLEC HECL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CLEC4C

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA855470LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA855470LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.